Institute of Medical & Biological Engineering

THE INFLUENCE OF ENVIRONMENTAL CONDITIONS ON THE WEAR OF AN ALL-POLYMER KNEE IMPLANT





Raelene M Cowie¹, Adam Briscoe², John Fisher¹, Louise M Jennings¹

¹Institute of Medical and Biological Engineering, School of Mechanical Engineering University of Leeds, UK ²Invibio Biomaterial Solutions Ltd., Lancashire, UK. Contact email: r.cowie@leeds.ac.uk



PEEK Optima®

Background

PEEK Optima® has been considered for use as an alternative arthroplasty bearing material due to its low wear rates, the low biological activity of its wear debris and clearance for clinical use [1].

This study investigated the potential to use PEEK Optima® as an alternative to cobalt chrome in the femoral component of a total knee replacement to give a metal free implant. The wear performance of all polyethylene tibial components was assessed under different environmental conditions.

Materials and Methods



Experimental wear simulation was carried out on 6 station ProSim knee simulators (Simulation Solutions, UK), one run at room temperature and the other at elevated temperature (~33°C).

- Test conditions used:
- 25% bovine serum in 0.03% sodium azide
- Kinematic conditions (Figures 2 & 3):
 - Axial Force (AF) up to ~2800N, Flexion extension (FE) 0-60°, Tibial Rotation (TR) ±5°, Anterior-Posterior Displacement (AP) 10mm (Leeds high kinematics) [3]
- . 5 million cycles (MC) at each temperature

Figure 1: Metal-free knee

- Six injection moulded PEEK Optima® femoral components (Invibio Biomaterial Solutions, UK), initial mean surface roughness (Ra) ~0.02µm (Figure 1)
- Six cobalt chrome femoral components, initial Ra ~0.02µm
- · Same design of femoral for both groups
- All polyethylene GUR1020 UHMWPE tibial components (conventional, ethylene oxide sterilised)
- . All components were cruciate retaining, right, mid-size
- Wear of the UHMWPE tibals was assessed gravimetrically
- Surface topography of the femoral components was measured using contact profilometry
- Lubricant temperature was measured using a thermocouple




Figure 2: Input AF and FE [2]





Figure 3: Input AP and TR [3]

Statistical analysis carried out using ANOVA with significance taken at p<0.05

Results

After 5MC of wear testing at room temperature, the mean wear rate of the conventional metal-on polyethylene implant was 2.6±1.5 mm³/MC and for the all-polymer knee was 4.2±5.4 mm³/MC, there was no significant difference in wear rate against the different materials p=0.27. At elevated temperature, the wear rates were lower against both materials at 0.2±0.5 and 1.5±2.5 mm³/MC respectively (Figure 4).



5 0 5 10 15 20 25 30 35 40 45

Cohalt Chrome

Figure 4: Mean wear rate ± 95% confidence limits of UHMWPE tibials (n=3)

Figure 5: Images of the femorals after 5MC of wear testing at room temperature (x10) taken using an Alicona G5

- At the conclusion of the studies, a high density of linear scratching was visible on the PEEK Optima® implants parallel to
 the principal direction of sliding (Figure 5). There was a significant increase (p<0.05) in the measured Ra of the PEEK
 femorals for both conditions; for the room temperature study, the post test Ra was 0.32±0.14µm.
- In the room temperature test, the mean bulk lubricant temperature was ~2°C higher in the all-polymer implant than the
 conventional materials

Discussion

- Both the conventional metal-on-polyethylene and all-polymer knee replacements produced low wear rates (<5mm³/MC).
- Testing at elevated temperature led to a reduction in wear rate of tibials against both femoral materials, likely due to protein
 precipitation and deposition on the surface of the implants [4].
- During wear testing, the PEEK Optima femorals became scratched however, this did not influence the wear rate which remained constant over the test duration.

Significance

- The wear performance of the all-polymer (PEEK-on-UHMWPE) knee was comparable to a conventional metal-onpolyethylene implant
- Environmental conditions such as lubricant temperature influence wear and tribology, testing at an elevated temperature introduced a test artefact which would not be seen in room temperature tests or in vivo.

Peteronee

[1] Kurtz, S.M. Biomaterials, 2007, 28(32):4845-69 [2] ISO 14243-3

[3] McEwen, H. M. J. J Biomechanics, 2005, 38(2):357-65[4] Wang, A. J Mater Sci Mater Med. 2012, 23: 1533-42

Financial Disclosur

JF is an NIHR Senior Investigator and his research is supported through the NIHR Leeds Musculoskeletal Biomedical Research Unit.

a kina włado omonie

This work was supported by Invibio Ltd and the Innovation and Knowledge Centre in Medical Technologies funded by the EPSRC, TSB and BBSRC. It was partially funded through WELMEC, a centre of Excellence in Medical Engineering funded by the Wellcome Trust and EPSRC under grant number WT 088908/Z/09/Z.

